English

## MR Conditional Disposable EEG [Electroencephalography] Cup Electrodes

### Intended Use

The MR Conditional Disposable EEG Cup Electrode is intended for use in the recording of the Electroencephalography [EEG], Evoked Potentials [EP] or as a Ground or Reference in an EEG or EP recording. This device is provided non-sterile for Single Patient Use Only and may remain on the patient in a MR environment under specific conditions.

### Sterilization

Non-sterile

### Caution

Federal [USA] law restricts this device to sale by or on the order of a physician and it should only be used in compliance with accepted industry standards. RhythmLink International, LLC is not responsible for injury, infection or other damage resulting from the use or misuse of this product.

MR Conditional Disposable EEG Cup Electrodes are for professional use only and should only be used in compliance with accepted industry standards. **The included extension cables [Fig. 1] are MR Unsafe.** Remove all extension cables before entering a MR environment.

### Packaging

Each Package Contains:

6, 12, 18 or 24 MR Conditional Disposable EEG Cup

Electrodes

6, 12, 18 or 24 MR Unsafe Extension cables

### Instructions for Use

Clean application site. Apply Cup Electrode using Weaver Ten20 conductive paste. MR Conditional Disposable EEG Cup Electrodes are only approved for use with Weaver Ten20 conductive paste. Remove all extension cables before entering an MR environment. When finished, remove electrodes and clean application sites.

### MRI Safety Information

Non-clinical testing has demonstrated that the MR Conditional Disposable EEG Cup Electrode [Fig. 2] is MR Conditional in configurations of 2 to 48 electrodes. These electrodes can safely remain on a patient during a MR scan for 15 minutes under the following conditions:

- Static magnetic field of 1.5 or 3.0 Tesla  
Maximum spatial gradient field of 4,000 gauss/cm [40T/m] or less
- Maximum whole-body averaged specific absorption rate of 2 W/kg in the Normal Operating Mode
- **Remove extension cables [Fig. 1] before entering an MR environment. They are MR Unsafe.**

Under the scan conditions defined above, the Disposable Cup EEG Electrode is expected to produce a maximum temperature rise of less than 1.4°C at 1.5 T and 3.1°C at 3 T after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the Cup Electrode extends less than 3 mm from the Cup Electrode when imaged with a gradient echo pulse sequence and a 3 T MRI system.

### RF Induced Heating Information

**1.5 Tesla Systems:** In non-clinical testing the MR Conditional Disposable EEG Cup Electrode, [in configurations of 2 to 48 electrodes] produced a temperature rise of less than 1.6°C [with a background temperature increase of ≈ 0.5°C] at a maximum head/whole body averaged specific absorption rate

[SAR] of ≈2.3 W/kg assessed by calorimetry for 15 min. of continuous MR scanning with whole body coil in a 64 MHz [1.5 Tesla equivalent] Medical Implant Test System, Zurich Medtech AG [Software: MITS-DUALBAND 1.2.5.2].

**3.0 Tesla Systems:** In non-clinical testing the MR Conditional Disposable EEG Cup Electrode, [in configurations of 2 to 48 electrodes] produced a temperature rise of less than 4.9°C [with a background temperature increase of ≈0.8°C] at a maximum head/whole body averaged specific absorption rate [SAR] of ≈3.2 W/kg assessed by calorimetry for 15 min. of continuous MR scanning with whole body coil in a 3 Tesla Magnetom Trio Siemens [Software: Numaris/4, syngo MR A30] MR Scanner.

**The MR Conditional Disposable EEG Cup Electrode, dual and multiple configuration [2 to 48 electrodes] have not been tested in simultaneous combination with other devices.**


### Artifact Information

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this device.

MR image artifacts can affect the device surrounding on each side from the device surface as follows:

Worst-case artifacts of	Spin Echo	Gradient Echo
Test object length	2.81 mm	2.52 mm
Test object diameter	2.28 mm	2.07mm

**The included extension cables are MR Unsafe.** Remove all extension cables before entering an MR environment.

[Fig. 1: MR Unsafe extension cable]

[Fig. 2: MR Conditional Disposable Cup EEG Electrode]

*This product has a useful shelf life of 5 years from the date of manufacture.*

RhythmLink International, LLC  
1140 First Street South  
Columbia, SC, USA 29209-3540  
+1.866.633.3754 [toll-free]  
+1.803.252.1222  
+1.803.252.1111 [fax]  
sales@rhythmink.com  
Rhythmink.com

**European Authorized Representative**  
Mdi Europa GmbH  
Langenhagener Str.71  
D-30855 Hannover-Langenhagen

**FOR SINGLE USE ONLY  
ASSEMBLED IN CHINA**

**Caution:** Federal Law [USA] restricts this device for sale by or on the order of a licensed medical practitioner. Rhythmink® is a registered trademark of RhythmLink International, LLC.

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Rhythmink® è un marchio commerciale registrato di Rhythmink International, LLC.